



SLIDE MISTER



1




Gently open the back cover by sliding it upwards.

4




Slide the cover downwards, the blue indicator will light up and mist will spray out of the spray nozzle.

2



Open tank rubber stopper.

3



Fill the tank with Copperine using the dropper

5



Aim towards the bottom of the nose while inhaling for 3-5 seconds. Suggested use: 3 times daily. May be used more frequently as needed, up to 10 times per day. Suitable for adults and children 1 year and older. Avoid spraying into the eyes.

Maintenance

- For best use and to increase the longevity of the product fully charge the Mister before using it.
- To keep your Mister clean, gently wipe the outside with dish soap and water. Do not rinse the Mister or soak it in water as it is not waterproof.
- Do not run the Mister without Copperine in the tank.
- If the mist is not spraying and the Mister has a adequate charge, give the Mister a very gentle shake and let it stand for a minute to clear any air bubbles.
- Keep the tank clean by removing it and cleaning it with distilled/ purified water - do not use soap

Copperine is for nasal use only and is not intended to replace a balanced diet or nutritional supplements. Not recommended for individuals with Wilson's Disease. These statements have not been evaluated by the Food and Drug Administration. **This product is not intended to diagnose, treat, cure, or prevent any disease.**